


ZOLL M2™

ZOLL®



**Empowering You to
Save More Lives**


More Trusted Technology, More Value

ZOLL defibrillators are designed to meet all your resuscitation needs. The ZOLL M2™ is no exception. From patient monitoring to resuscitation with high-quality CPR, the ZOLL M2 monitor/defibrillator embodies proven, reliable ZOLL core technology to ensure optimal patient care.



More Trusted Technology

As part of ZOLL's portfolio of resuscitation products, the ZOLL M2 monitor/defibrillator promotes consistent, high-quality CPR and high-current defibrillation for adults and children.



RECTILINEAR BIPHASIC DEFIBRILLATION


ZOLL's unique high-current waveform has been shown to be superior to monophasic waveforms for termination of ventricular and atrial fibrillation. It has been proven most effective in treating high-impedance patients when compared to truncated exponential waveforms at the same or higher energy levels.



PATENTED 40 MSEC PACING WAVEFORM

Unique, constant-current, 40 msec pacing has the highest capture rate at the lowest average current required, ensuring efficacy and patient comfort. Validated with over 4,000 patients in more than 16 studies.*

*Studies available upon request.



CPR FEEDBACK

Delivery of unaided CPR is not good enough. ZOLL has long recognized the importance of CPR feedback. Our Real CPR help® is used worldwide helping clinicians deliver high-quality CPR – the cornerstone of successful resuscitation.



REUSABLE HANDS-FREE SOFT PADDLES

ZOLL's trusted technology replaces external paddles and allows for shock delivery with Real CPR Help® included.

More Support

The ZOLL M2 professional defibrillator helps you deliver high-quality, resuscitation care to adult and pediatric patients by offering an array of life-saving features:

- CPR feedback to ensure delivery of high-quality CPR and successful resuscitation
- Cardioversion and pacing for patients with tachycardia and bradycardia
- AED mode with proven algorithms for both adults and children



More Value

Designed for clinicians, hospital staff, and emergency responders to offer high-quality patient care, the ZOLL M2 incorporates proven, best-in-class technology, including: defibrillation, pacing, cardioversion, and real-time CPR feedback—all for a lower total cost of ownership, empowering you to help save more lives.



More Scalable

From a simple shock box to the most complex patient management challenges, the ZOLL M2 offers various options for all patient resuscitation needs, including SpO₂, 12-lead ECG, NIBP, EtCO₂, Impedance Pneumography, and Temperature.

